



Dome® Controller Specifications

Make the *Right Choice*.

DISPLAY	Nvidia Graphics Board					
	NVS 290/295	FX380	Quadro 600	Quadro 2000D	Quadro 2000	Quadro 4000
GX2	XP, Win7	XP, Win7	XP, Win7	XP, Win7	XP, Win7 ³	XP, Win7 ³
Dome E2	XP, Win7	XP, Win7	XP, Win7	XP, Win7	XP, Win7 ³	XP, Win7 ³
Dome E2cHB	XP, Win7	XP, Win7	XP, Win7	XP, Win7	XP, Win7 ³	XP, Win7 ³
Dome E3	NS	XP, Win7	XP, Win7 ²	XP, Win7	XP, Win7 ³	XP, Win7 ³
Dome E3cHB	NS	XP, Win7	XP, Win7 ²	XP, Win7	XP, Win7 ²	XP, Win7 ²
Dome E3cHB/DP	XP, Win7 ¹	XP, Win7	XP, Win7	XP, Win7	XP, Win7	XP, Win7
Dome S3c	XP, Win7 ¹	XP, Win7	XP, Win7	XP, Win7	XP, Win7	XP, Win7
Dome E4c	NS	XP, Win7	XP, Win7 ²	XP, Win7	XP, Win7 ²	XP, Win7 ²
Dome E5	NS	XP, Win7	XP, Win7 ²	XP, Win7	XP, Win7 ³	XP, Win7 ³
Dome S10	NS	XP, Win7	XP, Win7 ²	XP, Win7	XP, Win7 ³	XP, Win7 ³

Nvidia Graphics Boards require driver version 275.89

All Graphics Boards support both 32 and 64 bit Operating Systems

Windows 7 support requires CXtra 6.2.3 or greater

¹ Supports Dual Head using DMS59 to Display Port Cable Adapter

² Supports Single Head only using the DVI connector, 2 boards are required for Dual Head

³ Display Port to DVI Adapter required for Display Port output(s)

NS = Not Supported



www.ndssi.com

Corporate Headquarters

5750 Hellyer Avenue
 San Jose, CA 95138 (USA)
 Tel: 408 776 0085
 Toll Free: 866 637 5237
 Email: info@ndssi.com

Europe

Nijverheidscentrum 28
 2761 JP Zevenhuizen (ZH)
 The Netherlands
 Tel: + 31 180 63 43 56
 Email: info-emea@ndssi.com

Asia Pacific

Takanawa Kaneo Bldg., 6F
 3-25-22 Takanawa
 Minato-ku Tokyo, Japan, 108-0074
 Tel: + 81 3 5475 1835
 Email: info@ndssi.jp