

ENDOVIEW™ 24

High-Definition Surgical Visualization



HD



 **EndoVue**

Fully compliant for medical use in surgery, the 24" EndoVue™ from NDSsi offers a high-quality, low-cost solution for surgical imaging environments. Its compact and lightweight design makes it very maneuverable and easy to install on a surgical cart or boom arm. The new 24" EndoVue™ is an affordable alternative to a full-featured surgical display but can still accommodate high-definition signals from a variety of medical imaging sources, including endoscopes, ultrasound, PACS, and vital signs. It also features a fanless cooling design to minimize the risk of spreading airborne contaminants near the sterile field.

- Purpose-Built for Use in the OR
- High-Quality, Low-cost Solution
- Lightweight, Highly-Maneuverable
- Programmable User Preferences
- Fanless Cooling Design

ENDOVIEW™ 24

High-Definition Surgical Visualization

HD



PRODUCT DETAILS

- Large, 24" Diagonal Size
- HD Imaging Compatible
- High Resolution (1920 x 1200)
- Programmable User and Modality Preferences
- Fanless Cooling, Sealed Front Enclosure
- Multiple Analog & Digital Video Inputs
- Safe Low Voltage 24 VDC Remote Power Supply
- Compact and Lightweight
- Flexible Standard VESA Mounting

SUPPORTED DEVICES/APPLICATIONS

SOURCE

DVI, HD-SDI
HD-RGBS, HD-YP_bP_r,
RGS, YP_bP_r, SDI
S-Video, Composite
VGA

APPLICATION

HD Endoscopy, PACS, Vital Signs, Room Camera
HD Endoscopy
SD Endoscopy
Fluoroscopy, SD Endoscopy, Ultrasound
PACS, Endoscopy, Ultrasound, Angiography

REGULATORY COMPLIANCE

UL 60601-1; EN 60601-1; EN 60601-1-2; MDD 93/42/EEC & 2007/47/EC; CE;
CAN/CSA C22.2 NO. 601.1; Class I Medical Device; FCC Part 15; RoHS; WEEE

NDSSI QUALITY SYSTEM COMPLIANCE

ISO 9001:2008 and ISO 13485:2003; FDA Registration # 2954921 & 1226517



Corporate Headquarters

5750 Hellyer Avenue
San Jose, CA 95138 (USA)
Tel: 408 776 0085
Toll Free: 866 637 5237
Email: info@ndssi.com

Europe

Nijverheidscentrum 28
2761 JP Zevenhuizen (ZH)
The Netherlands
Tel: + 31 180 63 43 56
Email: info@ndssi.nl

Asia Pacific

Takanawa Kaneo Bldg., 6F
3-25-22 Takanawa
Minato-ku Tokyo, Japan, 108-0074
Tel: + 81 3 5475 1835
Email: info@ndssi.jp