

NEW

EndoVue® 24"

Value Endoscopy Visualization



HD ENDOVIEW

- Optimized for endoscopy applications
- Outstanding image quality
- Multi-modality imaging capabilities
- Low power consumption and fanless cooling
- Robust and reliable
- Minimum cost of ownership

For the medical professional seeking to maximize their capital investment and upgrade their aging monitor, the EndoVue® series delivers the image performance of a high-end endoscopy display at a value price.

EndoVue is a versatile visualization platform, supporting multi-modality medical imaging requirements. It is plug-and-play compatible with endoscopic camera systems from leading manufacturers. Designed and built for medical use, it is FDA listed, TUV certified, and CE marked, offering patient safety, clinical performance, and durability compared to non-medical grade displays.

The EndoVue display accommodates standard analog and digital high-definition signals from a variety of medical imaging sources, including ultrasound, PACS, and patient monitor vital signs. It features a fanless cooling design to minimize the risk of spreading airborne contaminants within the sterile field, and industry-leading LED backlight technology for reduced power consumption.

EndoVue® 24"

Value Endoscopy Visualization

SPECIFICATIONS

Screen Size	24" diagonal
Resolution	1920 x 1200
Aspect Ratio	16:10
View Angle	178°
Luminance	300 cd/m ² typical
Color Gamut	BT.709 or SMPTE-C
Display Warranty	1 Year
Dimensions (WxHxD)	21.9" x 16.1" x 2.8"/566 x 408 x 70 mm
Weight	14.7 lbs/6.7 kg
Green Compliance	REACH, RoHS-2, WEEE
Regulatory Compliance	ANSI/AAMI ES60601-1, CAN/CSA C22.2 No. 60601-1, CAN ICES-3(B)/NMB-3(B), FCC Class B, EN60601-1, EN60601-1-2, CE, MDD 93/42/EEC, CCC, EAC, Conflict Minerals, Class I Medical Device

ORDER INFORMATION

EndoVue 24

90K0070	EndoVue 24" - North America
90K0071	EndoVue 24" - South America
90K0072	EndoVue 24" - Japan
90K0073	EndoVue 24" - China
90K0074	EndoVue 24" - Emerging Market
90K0075	EndoVue 24" - Europe
26B0102	EndoVue Display Stand - 19", 21", 24"



Global Headquarters

www.ndssi.com
5750 Hellyer Ave
San Jose, CA 95138
USA 🇺🇸
+1 408 776 0085
info@ndssi.com

Europe

GSI Group Europe GmbH
Münchner Strasse 2A
82152 Planegg
Germany 
+49 89 31 707 100
info.gsig-europe@gsig.com

Asia Pacific

East Square Omori
6-20-14 Minamioi, Shinagawa-ku
Tokyo 140-0013
Japan
+81 3 5753 2466
info@ndssi.jp

Oriental Media Center
Suite 2302, Tower C
No.4, Guang Hua Road
Chao Yang District
Beijing, 100026
China
+86 10 8559 7859
info@ndssi.com.cn

NDS Quality System

ISO 13485:2003
FDA Registration #2954921