

EndoVue® 24"

Value Endoscopy Visualization



ENDOVIEW

- Optimized for endoscopy applications
- Outstanding image quality
- Multi-modality imaging capabilities
- Low power consumption and fanless cooling
- Robust and reliable
- Minimum cost of ownership

For the medical professional seeking to maximize their capital investment and upgrade their aging monitor, the EndoVue® series delivers the image performance of a high-end endoscopy display at a value price.

EndoVue is a versatile visualization platform, supporting multi-modality medical imaging requirements. It is plug-and-play compatible with endoscopic camera systems from leading manufacturers. Designed and built for medical use, it is FDA listed, TUV certified, and CE marked, offering patient safety, clinical performance, and durability compared to non-medical grade displays.

The EndoVue display accommodates standard analog and digital high-definition signals from a variety of medical imaging sources, including ultrasound, PACS, and patient monitor vital signs. It features a fanless cooling design to minimize the risk of spreading airborne contaminants within the sterile field, and industry-leading LED backlight technology for reduced power consumption.

EndoVue® 24"

Value Endoscopy Visualization



SPECIFICATIONS

Screen Size	24" diagonal
Image Size	20.4" x 12.8"/518 x 324 mm (WxH)
Resolution	1920 x 1200
Aspect Ratio	16:10
View Angle	178°
Luminance	300 cd/m ² typical
Contrast Ratio	1,000:1 (typical)
Numbers of Colors	16.7 Million (8-bit)
Response Time	14ms
Latency	20ms
Video Inputs	DVI-D, RGBS/VGA/YpPr (HD-15 or DVI-I via adapter), S-Video, Composite
Video Output	DVI-D
Warranty	1 Year standard, optional extended 1 year
Dimensions	21.9" x 16.1" x 2.8"/566 x 408 x 70 mm (WxHxD)
Weight	14.7 lbs/6.7 kg
Shipping Weight*	22 lbs/9.9 kg
Shipping Dimensions	25" x 8" x 23"/64 x 21 x 59 cm
Green Compliance	REACH, RoHS-2, WEEE, Conflict Minerals
Regulatory Compliance	ANSI/AAMI ES60601-1, CAN/CSA C22.2 No. 60601-1, CAN ICES-3(B)/NMB-3(B), FCC Class B, EN60601-1, EN60601-1-2, CE, MDR 2017/745, CCC, EAC, Class I Medical Device
Supplied Accessories	Power Supply, 100-240 VAC, 50-60 Hz to 24 VDC, 2.7 A AC power cord, region appropriate

*Shipping weights are approximate.

ORDER INFORMATION

EndoVue 24

90K0070	EndoVue 24" - North America
90K0071	EndoVue 24" - South America
90K0072	EndoVue 24" - Japan
90K0073	EndoVue 24" - China
90K0074	EndoVue 24" - Emerging Market**
90K0075	EndoVue 24" - Europe

OPTIONAL ACCESSORIES

26B0102	Display Stand - 19", 21", 24"
35X0100	Fiber Converter 5 VDC Power Cable
73A1015	EndoVue 24" Extended 1 Year Warranty

**Emerging Market model eligible countries (non-eligible countries will void the warranty)

Bahrain | Bangladesh | Botswana | Bulgaria | Côte d'Ivoire | Croatia | Cyprus | Czech Republic | Egypt | Estonia | Ghana | Hungary | India | Indonesia | Jordan | Kenya | Lithuania | Macedonia | Malaysia | Malta | Mauritius | Morocco | Nigeria | Oman | Pakistan | Philippines | Poland | Qatar | Romania | Russia | Serbia | Slovakia | Slovenia | South Africa | Sri Lanka | Taiwan | Thailand | Tunisia | Turkey | UAE | Vietnam

Global Headquarters

www.ndssi.com
5750 Hellyer Ave
San Jose, CA 95138
USA 🇺🇸
+1 408 776 0085
info@ndssi.com

Europe

Novanta Europe GmbH
Parking 57-59
85748 Garching
Germany EC REP
+49 89 31 707 100
info.novanta-europe@novanta.com

Asia Pacific

Novanta Japan
East Square Omori
6-20-14 Minamioi, Shinagawa-ku
Tokyo 140-0013
Japan
+81 3 5753 2466
info@ndssi.com

Oriental Media Center

Suite 2302, Tower C
No.4, Guang Hua Road
Chao Yang District
Beijing, 100026
China
+86 10 8559 7859
info@ndssi.com

NDS Quality System

ISO 13485
FDA Registration #2954921