

# EndoVue® 32"

## Value Endoscopy Visualization



### ENDOVIEW

- Optimized for endoscopy applications
- Outstanding image quality
- Multi-modality imaging capabilities
- Low power consumption and fanless cooling
- Robust and reliable
- Minimum cost of ownership

For the medical professional seeking to maximize their capital investment and upgrade their aging monitor, the EndoVue® series delivers the image performance of a high-end endoscopy display at a value price.

EndoVue is a versatile visualization platform, supporting multi-modality medical imaging requirements. It is plug-and-play compatible with endoscopic camera systems from leading manufacturers. Designed and built for medical use, it is FDA listed, TUV certified, and CE marked, offering patient safety, clinical performance, and durability compared to non-medical grade displays.

The EndoVue display accommodates standard analog and digital high-definition signals from a variety of medical imaging sources, including ultrasound, PACS, and patient monitor vital signs. It features a fanless cooling design to minimize the risk of spreading airborne contaminants within the sterile field, and industry-leading LED backlight technology for reduced power consumption.

# EndoVue® 32"

## Value Endoscopy Visualization

### SPECIFICATIONS

Screen Size	31.5" diagonal
Image Size	27.5" x 15.5"/698 x 393 mm (WxH)
Resolution	1920 x 1080
Aspect Ratio	16:9
View Angle	178°
Luminance	400 cd/m <sup>2</sup> typical
Contrast Ratio	1,000:1 (typical)
Numbers of Colors	16.7 Million (8-bit)
Response Time	8ms
Latency	20ms
Video Inputs	DVI-D, RGBS/VGA/YpPr (HD-15 or DVI-I via adapter), S-Video, Composite
Video Output	DVI-D
Warranty	1 Year standard, optional extended 1 year
Dimensions	29.7" x 19.1" x 3.9"/753 x 485 x 98 mm (WxHxD)
Weight	21.3 lbs/9.7 kg
Shipping Weight*	29 lbs/13.1 kg
Shipping Dimensions	35" x 9" x 27"/89 x 23 x 69 cm
Green Compliance	REACH, RoHS-2, WEEE, Conflict Minerals
Regulatory Compliance	ANSI/AAMI ES60601-1, CAN/CSA C22.2 No. 60601-1, CAN ICES-3(B)/NMB-3(B), FCC Class B, EN60601-1, EN60601-1-2, CE, MDR 2017/745, EAC, CCC, Class I Medical Device
Supplied Accessories	Power Supply, 100-240 VAC, 50-60 Hz to 24 VDC, 2.7 A AC power cord, region appropriate

\*Shipping weights are approximate.

### ORDER INFORMATION

#### EndoVue 32

90K0080	EndoVue 32" - North America
90K0081	EndoVue 32" - South America
90K0082	EndoVue 32" - Japan
90K0083	EndoVue 32" - China
90K0084	EndoVue 32" - Emerging Market**
90K0085	EndoVue 32" - Europe

### OPTIONAL ACCESSORIES

35X0100	Fiber Converter 5 VDC Power Cable
73A1016	EndoVue 32" Extended 1 Year Warranty

\*\*Emerging Market model eligible countries (non-eligible countries will void the warranty)

Bahrain | Bangladesh | Botswana | Bulgaria | Côte d'Ivoire | Croatia | Cyprus | Czech Republic | Egypt | Estonia | Ghana | Hungary | India | Indonesia | Jordan | Kenya | Lithuania | Macedonia | Malaysia | Malta | Mauritius | Morocco | Nigeria | Oman | Pakistan | Philippines | Poland | Qatar | Romania | Russia | Serbia | Slovakia | Slovenia | South Africa | Sri Lanka | Taiwan | Thailand | Tunisia | Turkey | UAE | Vietnam



#### Global Headquarters

www.ndssi.com  
5750 Hellyer Ave  
San Jose, CA 95138  
USA 🇺🇸  
+1 408 776 0085  
info@ndssi.com

#### Europe

Novanta Europe GmbH  
Parking 57-59  
85748 Garching  
Germany   
+49 89 31 707 100  
info.novanta-europe@novanta.com

#### Asia Pacific

Novanta Japan  
East Square Omori  
6-20-14 Minamioi, Shinagawa-ku  
Tokyo 140-0013  
Japan  
+81 3 5753 2466  
info@ndssi.com

#### Oriental Media Center

Suite 2302, Tower C  
No.4, Guang Hua Road  
Chao Yang District  
Beijing, 100026  
China  
+86 10 8559 7859  
info@ndssi.com

#### NDS Quality System

ISO 13485  
FDA Registration #2954921